

AMMONIA STANDARD (NH₃ CAL)

LOT NUMBER: 855AM

SIZE: 1 x 5.5 ml

AMI015 GTIN: 05055273200256

AMI054 GTIN: 05055273200263

PREPARATION

Contents ready for use. Stable up to expiry quoted on the vial.

CONCENTRATION

0.05 g/l

304 µmol/l

0.52 mg/dl

SAFETY PRECAUTIONS

From Safety Data Sheet.

EC REP

Radox Teoranta, Meenmore,
Dungloe, Donegal,
F94 TV06, Ireland

21 Sep '23 me

THIS PAGE IS INTENTIONALLY BLANK